

BRM bromodomain TR-FRET Assay Kit

Cat. No. Kit-1935 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Size 9600 wells

Description

The SWI/SNF family of ATP-dependent chromatin remodeling enzymes catalyzes structural changes that allow proteins to access the nucleosomal DNA, alter the position of the nucleosomes on the DNA, or eject the histone octamer from the template. Mammalian SWI/SNF remodelers contain either BRM (SMARCA2) or BRG1 (SMARCA4) as their catalytic subunit. These subunits also contain a C-terminal bromodomain that binds to acetylated residues on histone tails. Small molecule inhibitors of BRM interactions with peptide binding partners hold promise as useful therapeutics for human disease. The BRM bromodomain TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format. This screening assay is robust, displaying a $Z'' > 0.6$.

Storage -80°C

Kit Components

BRM bromodomain Europium Chelate: 1 vial/420 wells; 5 vials/420 wells; 5 vials/2,100 wells; -80°C BRM bromodomain Ligand/APC Acceptor Mixture: 1 vial/420 wells; 5 vials/420 wells; 5 vials/2,100 wells; -80°C TR-FRET Assay Buffer (10X): 1 vial/2 ml; 1 vial/10 ml; 5 vials/10 ml; -20°C TR-FRET Assay Buffer Additive: 1 vial/200 mg; 1 vial/1 g; 5 vials/1 g; -20°C CH₃ Positive Control: 1 vial/2.5 nmol; 5 vials/2.5 nmol; 5 vials/12.5 nmol; -80°C 384-Well Solid Plate (low volume; black): 1 plate; 5 plates; 25 plates; RT Foil Plate Covers: 1 cover; 5 covers; 25 covers; RT

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA