

Phagocytosis Assay Kit (IgG FITC)

Cat. No. Kit-1996 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Size 1 ea

Description

The Phagocytosis Assay Kit (IgG FITC) employs latex beads coated with fluorescently-labeled rabbit IgG as a probe for the measurement of the phagocytic process in vitro. The engulfed fluorescent beads can be detected using a fluorescence microscope, allowing kinetic studies of phagocytosis at the single-cell level. In addition, the flow cytometric readout provides the advantage of visualizing perturbations in phagocytosis on the population level and, when combined with antibody staining, of specific cell types within complex populations. This kit provides enough Latex Beads-Rabbit IgG-FITC Complex for up to 750 samples.

Storage 4°C

Kit Components

Latex Beads-Rabbit IgG-FITC Complex: 1 vial/150 μ l 4°C Cell-Based Assay Buffer
Tablet: 2 tablets; RT Trypan Blue (10X): 1 vial/500 μ l RT

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA