

TAF1 bromodomains 1 and 2 TR-FRET Assay Kit

Cat. No. Kit-2033 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Size 9600 wells

Description

TAF1, a component of transcription factor IID, interacts with transcriptional activators to promote transcription. Recognition of acetylated lysine residues via its two bromodomains directs TAF1 binding to the core promoter sequences at the transcription start site. Through its interactions with the androgen receptor, TAF1 has been implicated in castration-resistant prostate cancer. TAF1 also has been shown to regulate expression of the pro-apoptotic protein p27kip1, potentially linking TAF1 expression to the cellular response to genotoxic stress. The TAF1 bromodomains 1 and 2 TR-FRET Assay Kit is a homogeneous TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format. This screening assay kit is robust, displaying a Z >0.6.

Storage -80°C

Kit Components

TAF1 bromodomains 1 and 2 Europium Chelate: 1 vial/420 wells; 5 vials/ 420 wells; 5 vials/ 2,100 wells; -80°C TAF1 bromodomains 1 and 2 Ligand/APC Acceptor Mixture: 1 vial/420 wells; 5 vials/ 420 wells; 5 vials/ 2,100 wells; -80°C TR-FRET Assay Buffer (10X): 1 vial/2 ml; 1 vial/10 ml; 5 vials/10 ml; -20°C TR-FRET Assay Buffer Additive: 1 vial/200 mg; 1 vial/1 g; 5 vials/1 g; -20°C H4 Positive Control: 1 vial/2.5 nmol; 5 vials/2.5 nmol; 5 vials/12.5 nmol; -80°C 384-Well Solid Plate (low volume; black): 1 plate; 5 plates; 25 plates; Room temperature Foil Plate Covers: 1 cover; 5 covers; 25 covers; Room temperature

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA