

Diamine Oxidase Assay Kit

Cat. No. Kit-2071 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The non-radioactive, fluorimetric DAO assay is based on the oxidation of putrescine to pyrroline plus NH ₃ and H ₂ O ₂ . The generated H ₂ O ₂ is then used by HRP to oxidize a dye making it fluorescent. The increase in fluorescence at $\lambda_{ex/em}$ = 530/585 nm is directly proportional to the enzyme activity.
Size	100 tests
Description	DIAMINE OXIDASE (DAO) also known as histaminase or amine oxidase (copper containing), is an enzyme involved in the metabolism, oxidation, and inactivation of histamine in animals. Highest content is observed in the digestive tract and placenta. An imbalance between histamine intake and the capacity for histamine degradation can lead to histamine intolerance (HIT). Measuring DAO activity in serum can be useful in diagnosing HIT.
Applications	For quantitative determination of diamine oxidase in biological samples (e.g. serum, plasma etc)
Storage	The kit is shipped at ambient temperature. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.
Kit Components	Assay Buffer: 10 mL HRP Enzyme: 120 μ Substrate: 120 μ H ₂ O ₂ Standard: 100 μ Dye Reagent: 120 μ
Detection method	FL530/585nm

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Compatible Sample Types	Serum, Plasma
Features & Benefits	• Fast and sensitive. Use of 10 µL sample. Linear detection range 0.5 to 6 U/L for 30 min reaction at 25°C. • Convenient. The procedure involves adding a single working reagent, and reading the fluorescence at 0 and 30 minutes. Room temperature assay. No 37°C heater is needed. • High-throughput. "Add-mix-read" type assay. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.
Assay time	30 min
Sensitivity	Detection Limit: 0.5 U/L

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA