# Recombinant Full Length Francisella Tularensis Subsp. Tularensis Nadh-Quinone Oxidoreductase Subunit K(Nuok) Protein, His-Tagged

**Cat. No.** RFL33819FF **Lot. No.** (See product label)

## SPECIFICATION

| | |
|---|---|
| **Product Overview** | Recombinant Full Length Francisella tularensis subsp. tularensis NADH-quinone oxidoreductase subunit K(nuoK) Protein (Q5NIM5) (1-110aa), fused to N-terminal His tag, was expressed in E. coli. |
| **Species** | Francisella tularensis |
| **Source** | E.coli |
| **ProteinLength** | Full Length (1-110) |
| **Form** | Lyophilized powder |
| **AA Sequence** | MRALKMNSISVSVTHGLIFSTLLFVISVAGIIINRRNILILLMSIELMLLAVNTNFLIFA NM HQQAMGGVFVFFIMAVAAAETAIGLAIVVAIFRKRKTIDLSKLNTLRG |
| **Purity** | Greater than 90% as determined by SDS-PAGE. |
| **Notes** | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| **Storage** | Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles. |

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA

| | |
|---|---|
| **Storage Buffer** | Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| **Reconstitution** | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20℃/-80℃. Our default final concentration of glycerol is 50%. Customers could use it as reference. |

## GENE INFORMATION

| | |
|---|---|
| **Gene Name** | nuoK |
| **Synonyms** | nuoK; FTT_0041; NADH-quinone oxidoreductase subunit K; NADH dehydrogenase I subunit K; NDH-1 subunit K |
| **UniProt ID** | Q5NIM5 |

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA