Amine Coupling Kit

Cat. No. Kit-0870 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	
Description	This Amine Coupling Kit is used with COOH or Amine Sensor Chips to covalently couple protein or antibody ligands to the surface. It has been optimized for effective and efficient coupling using a preconcentration technique which minimizes the amount of ligand needed to obtain sufficient immobilization levels. The blocking agent minimizes non-specific binding of the analyte. The Kit includes all necessary reagents for activating the COOH groups on the sensor surface (COOH Sensor Chips) or on the protein surface (Amine Sensor Chips) with EDC/NHS, immobilizing the protein ligand, and blocking the surface.
Size	10 immobilizations
Kit Components	EDC: 1 tube NHS: 1 tube Activation Buffer Blocking Solution

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA