Creatine Kinase (CK) Activity Colorimetric Assay Kit

Cat. No. Kit-0230 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The kit can detect creatine kinase activity less than 1 mU. (100 assays) Simple procedure; takes ~ 30-40 minutes;Sensitive and high-throughput adaptable;Kit contains all necessary reagents for measuring creatine kinase activity.
Size	100 assays
Description	Creatine Kinase (CK) also known as creatine phosphokinase (CPK) and ATP: creatine N-phosphotransferase is a common cellular enzyme (EC 2.7.3.2). It catalyzes the reversible conversion of creatine and ATP into ADP and phosphocreatine. CK is widely expressed in various tissues and cell types, with highest activity in striated muscles, heart tissue and brain. CK consists of two subunits: M (muscle) and B (brain), and has three isoenzymes: CK-MM (skeleton muscle), CK-MB (cardiac muscle), and CK-BB (brain). Increased CK level is associated with many diseases such as myocardial infarction, muscular dystrophy, pulmonary infarction and brain tumors. Accurate measurement of CK is crucial for early diagnosis, prediction and therapeutic strategy. In Creatine Kinase Activity Colorimetric Assay kit, creatine kinase converts creatine into phosphocreatine and ADP. The generated phoshocreatine and ADP reacts with CK Enzyme Mix to form an intermediate, which reduces a colorless Probe to a colored product with strong absorbance at 450 nm. The CK Activity Assay is high-throughput adaptable, simple and sensitive. This assay kit can detect Creatine Kinase activity less than 1 mU.
Applications	The kit can detect creatine kinase activity less than 1 mU.

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA

Usage	This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.
Kit Components	CK Assay BufferCK SubstrateATP (Lyophilized)CK Enzyme Mix (Lyophilized)CK Developer (Lyophilized)NADH Standard (Lyophilized)Positive Control (Lyophilized)
Detection method	Colorimetric
Compatible Sample Types	Animal TissueCell Culture SupernatantsPlasmaSerum

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA