Maleimide Green Fluorescence Assay Kit

Cat. No. Kit-0534 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Maleimide Assay Kit (Green Fluorescence) is a one-step fluorometric method for the quantification of maleimide groups in a variety of molecules such as proteins.
Description	Sensitive assays of maleimide and thiol groups are required for the efficient conjugation of proteins that are expensive and available only in small amounts. A variety of crosslinking reagents with a maleimide group are widely used for crosslinking proteins to proteins or proteins to other biomolecules. There are few reagents or assay kits available for quantifying the number of maleimide groups that are introduced into the first protein. All the commercial kits have tedious protocols.
Applications	Maleimide Assay Kit (Green Fluorescence) uses a proprietary dye that has enhanced fluorescence upon reacting with a maleimide.
Usage	For research use only (RUO)
Storage	Store the kit at -20°C and avoid from light.
Kit Components	Component A: Maleimide Green 1 vialComponent B: Reaction Buffer 500 µLComponent C: Assay Buffer 25 mLComponent D: N-ethylmaleimide Standard (10 mM) 50 µLComponent E: DMSO 200 µL
Detection method	Fluorometric
Features & Benefits	Broad Application: Can be used for quantifying maleimide group in a variety of molecules such as proteins.Sensitivity: Detects as low as 10 picomoles of

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA

	maleimide.Continuous: Easily adapted to automation without a separation step.Convenient: Formulated to have minimal hands-on time. No wash is required.Non-Radioactive: No special requirements for waste treatment.

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA