Oxalate Decarboxylase Activity Colorimetric Assay Kit

Cat. No. Kit-0641 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Oxalate Decarboxylase Activity Assay Kit (Colorimetric) is used to measure Oxalate Decarboxylase activity.
Description	Oxalate Decarboxylase (OXDC, EC 4.1.1.2) belongs to the cupin superfamily and is composed of two β-barrel domains. OXDC catalyzes the conversion of Oxalate into Formate and CO2 and plays an important role in stress response. In humans, high levels of oxalate can lead to various health problems including hyperoxaluria, kidney stones, and renal failure. Wood rotting fungi generates high levels of oxalate, causing rot in many crops including lettuce, soybean, dry bean, and tomato etc. Recent studies show that overexpression of OXDC in plants (e.g. tomato, soybean, lettuce, and tobacco) results in transgenic plants with resistance to fungal pathogenesis. Accurate measurement of Oxalate Decarboxylase activity is useful for a variety of therapeutic, diagnostic, and mechanistic studies. The Oxalate Decarboxylase Activity Assay Kit (Colorimetric) provides a quick and easy way to measure Oxalate Decarboxylase activity in various samples. In this assay, Oxalate Decarboxylase converts oxalate to formate, which subsequently converts a nearly colorless probe to a colored product with strong absorbance at 450 nm. The assay is simple, sensitive, high-throughput adaptable and can detect less than 20 μU of oxalate decarboxylase activity in a variety of samples.
Applications	Measurement of Oxalate Decarboxylase activity in various plant tissues.Analysis of cell signaling pathways such as glyoxylate and dicarboxylate metabolism
Target Species	Fungi, Plants

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA

Usage	For research use only (RUO)
Storage	Store the kit at -20°C, protected from light.
Kit Components	OXDC Assay Buffer I 20 mLOXDC Assay Buffer II 15 mLOXDC Substrate (Lyophilized) 1 vialOXDC Enzyme Mix (Lyophilized) 1 vialOXDC Probe (Lyophilized) 1 vialOXDC Positive Control 0.1 mLFormate Standard (100 mM) 0.1 mL
Detection method	Colorimetric
Compatible Sample Types	Fungi, Plant tissues

Tel: 1-631-559-9269 1-516-512-3133
Email: info@creative-biomart.com Fax: 1-631-938-8127
45-1 Ramsey Road, Shirley, NY 11967, USA